
TECHNICAL DATA SHEET

MI-KV-KIT

MI KNEE OPERATED VALVE KIT







DESCRIPTION

The MacDonald Industries Knee Operated Valve is another part of the MDI hands free tapware range. These valves are extremely robust and long lasting. These are a timed flow valve with a 7 second runtime but can also be held in the ON position to give the user a greater time of use.

Hygienic hands free activation for use in hospitals, kitchens and laboratories

- 15mm BSP inlet and outlet
- 17 second run time
- 150-350kPa water inlet pressure
- 0-60 degree water temperature

MODELS

MI-KV-KIT – Knee operated valve kit, 7 second run time

KIT INCLUDES:

- 1 x 150mm Gooseneck outlet
- 2 x 15mm Mini Stops
- 2 x 15mm BSP raised cover plates
- 2 x 15mm C/P Kinco Nuts
- 2 x 15mm Nylon Olives
- 1 x 15BSP C/P Copper Sink Tee
- 1 x 300 FxF Flexi Hose
- 1 x 450 FxF Flexi Hose



Melbourne (Head Office)

4/168 Jersey Drive, Epping Victoria 3076 sales@macdonaldindustries.com.au macdonaldindustries.com.au